

DRY / STORAGE BOX



PRODUCT PROTECTION

These boxes are designed to store samples and protect them against dust or neutral atmosphere (gas scavenging).

■ DIMENSIONS (mm)

TYPE	HEIGHT (H)	WIDTH(W)	DEPTH(D)
B1	300	300	300
B2	500	300	300
B3	750	300	300

■ CHARACTERISTICS

- Methyl methacrylate structure
- Door with neoprene seal
- Quarter turn closure
- 1 shelf by compartment

■ OPTIONS

- Laminated wood structure
- Inactinic structure, antistatic
- Gas passage, tap
- Flowmeter

■ FLAT-PAST VERSION

SAS pass-throughs allow the introduction of equipment directly into the clean room.

Standard dimension : 600x600x600 mm
Options : condamnation des portes (automatisme)

Note: for any specific requests ,consulte with our team of specialist



Gas sweeping cabinet



Sas with a lock

COMMISSIONING AND MAINTENANCE PROVIDED BY ADS LAMINAIRE

INF - COM - 13 - Rév1 - "Document non contractuel. Nous nous réservons la possibilité de faire évoluer les caractéristiques de nos produits sans préavis"