



ADS Laminaire

Air Dépoussiéré Stérile Laminaire

Product sheet



ADJUSTABLE ELECTRIC BASE



Description :

The 2-motor electric system offers exceptional flexibility and ergonomics.

Height adjustable on 2 adjustment motors, “ T-motion ”, it guarantees smooth and silent movement. The box controls 3 memories to save your favorite configurations.

The collision avoidance system “ T-touch ” ensures your safety.

Material(s):

High quality steel for remarkable strength and durability.

The top can be made to measure and available in several materials: PVC, TRESPA, INOX, PPH, etc.

Colors:


White RAL9010, Gray RAL9006, Black RAL9005

Specificities :

Adjustable height: From 625 mm to 1275 mm

Maximum load: 120 kg (Minus the weight of the tray, i.e. 80n to 110 kg depending on the tray). kind

Base plans (seen from above)



Base perfectly suited to our ranges

OPTIGEL
Cpure
IVGEL/IVGELSPACE
MINISPACE/MINIPURE
BIOII



**C-PURE + presentation video
Electric base**