

Complies with
ISO 14 644



MINISPACE

Vertical laminar flow cabinet, ISO 5 class

Protection of the product

Vertical ISO 5 sterile airflow

Fully automatic

ECM regulation and LED lighting

Small footprint

800 mm overall

Ultra-Cleanliness Solutions

Protect the product · Protect the operator · Protect the environment

 French manufacturer since 1968 · adslaminaire.com

PURPOSE

The **MINISPACE** consists of a supply and filtration unit that guarantees sterile working conditions. **ISO 5** class hoods protect the product by means of a vertical sterile airflow.

Fully automatic and fitted with LED lighting, it suits any laboratory layout — including short benches where a standard model will not fit.

KEY FEATURES

- Automatic ECM regulation
- Airflow velocity display
- Fully open front
- Transparent work chamber
- Stainless steel work surface
- LED lighting
- Ergonomics — side openings and inclined front



A format that changes everything

At 800 mm overall width and 510 mm deep, the MINISPACE goes where bench space is tight, without giving up ISO 5.

■ OPERATING PRINCIPLE

The supply and filtration unit delivers a **vertical** sterile airflow across the whole work surface. The air velocity and the positive pressure protect the product: ambient particles cannot enter the protected zone.

ECM regulation compensates for the gradual loading of the HEPA filter and keeps the airflow constant throughout its life.

■ DESIGN DETAILS

Side access

The side openings and the inclined front keep the view and the access clear, which offsets the cabinet's small format.

Velocity display

The airflow velocity is displayed at all times, so a glance is enough to confirm that the cabinet is working within specification.

DIMENSIONS AND SPECIFICATIONS

Dimensions (mm)	MINISPACE
	internal / overall
Width	780 / 800
Height	650 / 935 — on stand: 1 766
Depth	500 / 510
Weight	50 kg
Regulation	Automatic ECM, constant airflow
Lighting	LED

Source: ADS Laminaire general catalogue 2025-2026.

Fiche MINISPACE — Not a contractual document. We reserve the right to modify the specifications of our products without notice.

DOCUMENTS SUPPLIED ON DELIVERY

- Certificate of conformity
- Warranty certificate
- User manual in English

Commissioning and maintenance are carried out by ADS Laminaire, the manufacturer, throughout France and Europe.

OUR COMMITMENT

Engineering, manufacturing, test laboratory and after-sales service, all under one roof at our Aulnay-sous-Bois site.

A company certified **ISO 9001:2015** since 2008 and awarded the **NF 095** label since 1997. The first French company to specialise in laminar flow, since 1968.