ADS Laminaire
air dépoussiéré stérile laminaire
FABRIQUÉ EN FRANCE

ISO 9001:2015 CERTIFIED COMPANY

MADE IN FRANCE

Compliant with ISO 14 644



MINISPACE

Vertical laminar flow cabinet, ISO class 5

Product protection
Sterile vertical ISO 5 airflow

Fully automatic
ECM control and LED lighting

Small footprint
800 mm overall width

Ultra-Clean Solutions

Protect the product · Protect the user · Protect the environment

French manufacturer since 1968 · adslaminaire.com

INTENDED USE

The **MINISPACE** consists of a blower and filtration unit that ensures work in a sterile environment. **ISO 5** cabinets are designed to protect the product with a vertical sterile airflow.

Fully automatic and fitted with LED lighting, it suits any laboratory layout — including short benches where a standard model will not fit.



FEATURES

- **Automatic ECM control**
- **Airflow velocity display**
- **Plexiglas sash** — on the front
- **Transparent work area**
- **Stainless steel work surface**
- **LED lighting**
- **Ergonomics** — open sides and tilted front

A format that changes everything

With an overall width of 800 mm and a depth of 510 mm, the MINISPACE fits where bench space is limited, without giving up ISO class 5.

Ultra-Clean Solutions

Designed and made in France

34 rue Maurice de Broglie — ZI Les Mardelles
93600 Aulnay-sous-Bois, France
Tel. +33 (0)1 48 17 85 10 · Fax +33 (0)1 48 17 85 11
info@adslaminaire.com · adslaminaire.com

OPERATING PRINCIPLE

The blower and filtration unit delivers a sterile **vertical** airflow over the entire work surface. The airflow velocity and positive pressure protect the product: ambient particles cannot enter the protected zone.

ECM control compensates for gradual HEPA filter clogging and keeps the airflow constant throughout the filter's life.

SPECIAL FEATURES

Side ergonomics

The open sides and tilted front improve visibility and access, making up for the compact size of the cabinet.

Velocity display

The airflow velocity is displayed at all times, so you can see at a glance that the cabinet is working within its specifications.

Ultra-Clean Solutions

Designed and made in France

34 rue Maurice de Broglie — ZI Les Mardelles
93600 Aulnay-sous-Bois, France
Tel. +33 (0)1 48 17 85 10 · Fax +33 (0)1 48 17 85 11
info@adslaminaire.com · adslaminaire.com

DIMENSIONS AND SPECIFICATIONS

Dimensions (mm)	MINISPACE
	internal / overall
Width	780 / 800
Height	650 / 935 — on stand: 1 766
Depth	500 / 510
Weight	50 kg
Control	Automatic ECM, constant airflow
Lighting	LED

Source: ADS Laminaire general catalogue 2025-2026.
MINISPACE data sheet — Non-contractual document. We reserve the right to change the specifications of our products without notice.

DOCUMENTS SUPPLIED ON DELIVERY

- Certificate of conformity
- Warranty certificate
- User manual in French

Commissioning and maintenance carried out by the manufacturer, ADS Laminaire, throughout France and Europe.

OUR COMMITMENT

Design, manufacturing, test laboratory and after-sales service all on our Aulnay-sous-Bois site.

Company certified **ISO 9001:2015** since 2008 and **NF** labelled since 1997. A French company specialising in laminar flow since 1968.

Ultra-Clean Solutions

Designed and made in France

34 rue Maurice de Broglie — ZI Les Mardelles
93600 Aulnay-sous-Bois, France
Tel. +33 (0)1 48 17 85 10 · Fax +33 (0)1 48 17 85 11
info@adslaminaire.com · adslaminaire.com